

Sterile, Disposable Safe-Toe[®] Podiatry Burs

Featuring Rounded Bevels to Protect Patients' Skin

Safe-Toe[®] Podiatry Burs are specifically designed to prevent injuries to the skin and nail bed making them perfect for the management of at-risk patients.

TECHNICAL SPECIFICATIONS

- Single-use to guarantee a clean, sharp bur every time
- Surgical grade stainless steel shaft
- Consistent 2.35 mm shaft diameter
- Made in Europe
- Sterilized with gamma irradiation
- All natural diamond head
- Sold in quantities of 40 EA/CA

TOP GRINDER

111-MD30S

Selective treatment of diseased, thickened or deformed nails:

- Diamond coating on top only
- Smooth, non-coated side



SIDE GRINDER

111-MD31S

Ideal for a safe therapy at nail borders and close to nail beds:

- Diamond coating on side only
- Smooth, non-coated top



ROUND-EDGE GRINDER

111-MD32S

Perfect choice for skin removal, grinding and smoothing:

- Diamond coating on top and side
- Perfect for scarfskin and nail surfaces



RECOMMENDED ROTATION SPEED: 5,000-10,000 RPM (MIN-1)

We recommend a decrease in rotation speed during dry use and on at-risk patients in order to prevent increased temperature.