

steriSTAT® STERILE DISPOSABLE PODIATRY BURS

TAA Compliant • FDA Approved • ISO Certified Single-Use • Individually Packaged and Labeled
Compatible with Most Drills • Stainless Steel Shank • Standard Shaft Size • Sterilized with Gamma Irradiation

NATURAL DIAMOND

(Sold in Quantities of 40 EA / CA)

111-MD23S



Bud Shaped

Working L = 11.5 mm
Head Ø = 7.3 mm

111-MD24S



Pear Shaped

Working L = 12.1 mm
Head Ø = 7.4 mm


111-MD25S



Small Taper

Working L = 10.5 mm
Head Ø = 2.1 mm

111-MD26S



Small Short Barrel

Working L = 6.5 mm
Head Ø = 3.1 mm

111-MD27S



Olive Shaped

Working L = 13.8 mm
Head Ø = 7.4 mm

111-MD28S



Large Barrel

Working L = 11.6 mm
Head Ø = 7.4 mm

111-MDC2BS

(Qty. 25 EA / CA)



Callus Cutter

Working L = 15.7 mm
Head Ø = 8.2 mm

CARBIDE

(Sold in Quantities of 25 EA / CA)

111-MDCXCCS



Extra Cutter, Coarse

Working L = 14.4 mm
Head Ø = 6 mm

111-MDC2XCCS



Cutter, 2x Coarse

Working L = 14.5 mm
Head Ø = 5 mm

111-MDCC6S



Pear Shaped, Coarse

Working L = 10 mm
Head Ø = 5.7 mm

111-MDC7S



Small Cone, Fine

Working L = 12.9 mm
Head Ø = 6 mm

111-MDCC7S



Small Cone, Coarse

Working L = 12.9 mm
Head Ø = 6 mm

111-MDCC8S



Medium Cone, Coarse

Working L = 12.5 mm
Head Ø = 5.8 mm

111-MDCC9S



Large Cone, Coarse

Working L = 14.7 mm
Head Ø = 6 mm